

easi-lift shower stretcher



Easy adjustment

The Easi-Lift Shower Stretcher has been designed to help you move the user up and down with ease whilst they are on the stretcher. It is electronically operated with a height range of 500mm, which allows a safe and suitable working height to prevent the risk of back injury to the carer.

Contemporary design

The unobtrusive and neutral design of the Easi-Lift Shower Stretcher allows it to fit into a wide range of environments. The need for floor fixing has been eliminated making it the perfect choice for areas where the floor type, surface pipe work or cabling prevent floor installation.

Easy storage

The compact design coupled with the ability to fold the stretcher away when it is not in use makes this the ideal product for times when space is limited.

easi-lift shower stretcher

key features



- Electrically operated
- No floor fixing required
- Ideal for a wet room environment
- Stainless steel frame with powder coated finish
- Available in 6 sizes
- Length of stretcher available in 100mm increments from 1400mm – 1900mm
- Weight capacity of 200kg
- Travel range of 500mm
- Folds away to save on space
- Easy to clean surfaces
- Optional integral vinyl headrest
- Optional 1400mm safety guard

Stainless Steel Powder Coated		Size	User Weight kg	Table Dimensions			Travel Range H mm	Wall Frame Width A mm	Wall Frame Depth B mm	Safety Guard Height Y mm	Stored Width S mm
Without Headrest	With Headrest			L mm	W mm	M mm					
3191	3291	1	200	1400	650	1000	450 - 950*	750	110	145**	270
3192	3292	2	200	1500	650	1000	450 - 950*	750	110	145**	270
3193	3293	3	200	1600	650	1000	450 - 950*	750	110	145**	270
3194	3294	4	200	1700	650	1000	450 - 950*	750	110	145**	270
3195	3295	5	200	1800	650	1000	450 - 950*	750	110	145**	270
3196	3296	6	200	1900	650	1000	450 - 950*	750	110	145**	270

*dependant on installed height **safety guard is optional accessory

